

The logo for REZEX, featuring the word "REZEX" in a bold, sans-serif font. The letters "R", "E", and "Z" are green, while "E", "X", and "E" are purple. A small "TM" trademark symbol is located at the top right of the "X".

REZEX™

Reproducible Separation of Carbohydrate, Oligosaccharide, and Organic Acid Analysis

- Accurate and reproducible analysis
- Long column lifetimes



The Phenomenex logo, consisting of a stylized white circle with a registered trademark symbol (®) above it.

phenomenex®
...breaking with traditionSM



Trust Rezex™ For Your HPLC Analysis Of Carbohydrates, Oligosaccharides, And Organic Acids

Phenomenex Rezex HPLC ion-exclusion and ion-exchange columns are guaranteed to give you the performance you need. From drug formulation analysis to quality control testing of finished food products, Rezex columns consistently provide accurate and reproducible results. Try Rezex columns today risk-free and be rewarded with long lasting quality.

Try Rezex Risk Free! Rezex is a guaranteed alternative to:

- Bio-Rad® Aminex®
- Waters® Sugar-Pak™
- Supelco® SUPELCOGEL™
- Transgenomic® CARBOSep™

guarantee

If you are not completely satisfied with the performance of any Rezex column, as compared to a competing product of the same size and phase, simply return the Rezex column with your comparative data within 45 days for a FULL REFUND.



For Excellent Resolution And Reproducibility Of Sugars, Starches, Organic Acids, And Oligonucleotides

Broad Range of Phases to Perfectly Suit Your Application Needs	2
Column Selection - Variety Gives You the Power of Optimization	3
See the Difference! The Rezex Performance Advantage	
Sharper Peak Shape = Easy & Accurate Quantitation	4
Lower Backpressure = Longer Column Lifetimes & Faster Run Time Capability	5
Higher Efficiency = Baseline Separation of Critical Sample Components	5
The Proof is in the Performance – Select Applications	
Food and Beverages.....	6
Fermentation / Bioethanol.....	6
Carbohydrates / Saccharides.....	7-8
Alcohols.....	8
Other Applications	9
Specifications and Operating Recommendations	10
Ordering Information	11



Broad Range Of Phases To Perfectly Suit Your Application Needs

Rezex™ columns achieve separations based on multiple modes of interaction including ion-exchange, ion-exclusion, size exclusion, reversed phase, and partition. In many cases, this allows for separation of multiple compound classes with one column. The 8 % cross-linked materials exhibit excellent strength and durability while maintaining low system backpressures not typical with highly cross-linked resins. The 4 % cross-linked silver and sodium phases have a larger effective pore volume and are ideal for larger starches, such as oligosaccharides.



Table 1

			Cross Referencing Guide			
Phenomenex Rezex Phase Available	Description	Applications	Bio-Rad® Aminex®	Supelco® SUPELCOGEL™	Waters® Sugar-Pak™	Transgenomic®
RCM-Monosaccharide (L19 packing)*	8 % cross-linked CALCIUM	<ul style="list-style-type: none"> Monosaccharides and sugar alcohols, including sorbitol and mannitol from sweeteners and corn and cane sugars Class separation of di-, tri-, and tetra- saccharides 	HPX-87C 125-0095	SUPELCOGEL Ca	Sugar-Pak 1	CARBOsep™ CHO-820
RHM-Monosaccharide (L17 packing)*	8 % cross-linked HYDROGEN	<ul style="list-style-type: none"> Monosaccharides in combination with organic acids, fatty acids, alcohols, ketones, neutral compounds, or inorganic salts 	HPX-87H 125-0140	SUPELCOGEL C-610H & H	N/A	ICSep ION-300
RAM-Carbohydrate	8 % cross-linked SILVER	<ul style="list-style-type: none"> Unique selectivity complimentary to other Rezex column types 	N/A	N/A	N/A	–
RSO-Oligosaccharide	4 % cross-linked SILVER	<ul style="list-style-type: none"> High resolution of oligosaccharides up to 18 degrees of polymerization (Dp) 	HPX-42A 125-0097	SUPELCOGEL Ag1 & Ag2	N/A	–
RNO-Oligosaccharide	4 % cross-linked SODIUM	<ul style="list-style-type: none"> High resolution of oligosaccharides 	N/A	N/A	N/A	CARBOsep COREGEL-87N
RPM-Monosaccharide (L34 packing)*	8 % cross-linked LEAD	<ul style="list-style-type: none"> 100 x 7.8 mm dimensions available for sugar alcohol analysis according to the USP procedure Monosaccharides and sugar alcohol analysis Cellobiose, glucose, xylose, arabinose, and mannose and other cellulose products 	HPX-87P 125-0098	SUPELCOGEL Pb	N/A	CARBOsep COREGEL-87P
RNM-Carbohydrate (L58 packing)*	8 % cross-linked SODIUM	<ul style="list-style-type: none"> Analysis of mono-, di-, and tri- saccharides For matrices which contain high concentration of inorganic sodium (i.e.-molasses) 	HPX-87N 125-0143	N/A	N/A	–
ROA-Organic Acid (L22 packing)*	8 % cross-linked HYDROGEN	<ul style="list-style-type: none"> Organic acids Organic acids in combination with carbohydrates, alcohols, fatty acids, or neutral compounds Amino sugars 	HPX-87H 125-0140	SUPELCOGEL C-610H & H	N/A	–
RFQ-Fast Acid	8 % cross-linked HYDROGEN	<ul style="list-style-type: none"> Rapid screening of fruit quality Ethanol, acetic acid, glycerol, & standard alcohol mixtures 	Fast Acid 125-0100	N/A	N/A	–
RKP-Potassium	8 % cross-linked POTASSIUM	<ul style="list-style-type: none"> Analysis of glyphosate 	HPX-87K 125-0142	SUPELCOGEL K	N/A	CARBOsep COREGEL-87K
RCU-USP Sugar Alcohols (L19 packing)*	8 % cross-linked CALCIUM	<ul style="list-style-type: none"> Sugar analysis according to the USP procedures on a 250 x 4.0 mm dimension 	Sugar Alcohols 125-0094	N/A	N/A	–

* United States Pharmacopeia (USP)

Column Selection – Variety Gives You The Power Of Optimization

Rezex™ columns utilize different ligands. Selecting the right column will allow you to achieve the best resolution for your peaks of interest.

Column Selection Guideline:

1. Application examples (Table 1)
2. Cross referencing guide (Table 1)
3. Retention time chart (Table 2)
4. Our support staff – Knowledgeable consultants who can provide excellent recommendations on applications



Retention Times For Some Carbohydrates And Sugar Alcohols

Still not sure what phase to choose? It is important to ensure that you get adequate resolution between your peaks of interest. The retention time chart below gives the retention times for common carbohydrates and sugar alcohols on the most popular Rezex phases.

Table 2

Counter Ion Analyte	RAM Ag ⁺	RCM Ca ⁺²	RNM Na ⁺	RHM H ⁺	RPM Pb ⁺²
Adonitol (Ribitol)	11.54	14.93	11.10	11.11	20.15
D-Altrose	11.95	12.71	11.45	10.21	15.82
D-(-)-Arabinose	13.01	13.56	12.65	11.24	16.47
D-(+)-Cellobiose	8.86	8.60	8.49	8.02	11.00
D-(+)-Digitoxose	11.90	13.82	11.39	12.59	15.32
Dulcitol	11.64	21.61	11.10	10.71	33.25
Meso-Erythritol	12.31	15.49	11.78	12.14	19.82
D-(-) Fructose	12.05	13.65	11.76	10.31	17.71
L-(-)-Fucose	12.75	13.19	12.30	11.65	16.19
D-(+)-Galactose	11.87	11.73	11.47	10.19	14.94
Gentiobiose	8.70	8.40	8.40	7.87	10.53
D-(+)-Glucose	11.04	10.37	10.71	9.62	12.92
Inositol	12.59	13.35	12.14	9.98	18.87
Isomaltose	9.11	8.74	8.76	8.02	11.28
Lactose	9.27	9.03	8.78	8.32	11.89
Lactulose	9.75	10.32	9.23	8.57	13.95
D- Lyxose	12.41	14.06	11.98	10.68	16.66
D- Maltose	9.16	8.81	8.75	8.18	11.59
Maltotriose	8.27	8.10	7.94	7.51	11.02
Maltulose	9.25	9.47	8.82	8.27	12.40

* Partial hydrolysis results.

Counter Ion Analyte	RAM Ag ⁺	RCM Ca ⁺²	RNM Na ⁺	RHM H ⁺	RPM Pb ⁺²
D- Mannitol	11.36	17.82	10.80	10.59	24.90
D-(+)-Mannose	12.04	12.04	11.54	10.16	16.39
Melibiose	9.26	9.04	8.82	8.14	11.97
D-(+)-Melezitose	8.00	7.93	7.66	7.54*	9.94
D-(+)-Raffinose	8.10	8.16	7.76	7.88*	10.28
L-(+)-Rhamnose	11.50	12.18	11.00	10.90	14.47
D-(-)-Ribose	14.59	23.38	14.34	11.42	33.48
Salicin	18.51	18.58	17.36	14.98	26.81
D-Sorbitol	11.91	22.45	11.39	10.83	35.97
Stachyose	7.60	7.59	7.30	7.27	9.72
Sucrose	9.03	8.71	8.65	9.24*	11.00
Trehalose	8.91	8.72	8.49	8.32	11.01
Xylitol	12.69	22.01	12.16	11.78	32.38
D-(+)-Xylose	12.06	11.62	11.68	10.24	13.84


Conditions

Dimensions: 300 x 7.8 mm
Mobile Phase: Water (degassed)
Flow Rate: 0.6 mL/min
Temperature: 80 °C
Detection: RI @ 40 °C

Phenomenex® Rezex™ vs. BIO-RAD® AMINEX®


Sharper Peak Shape = Easy & Accurate Quantitation

► Saccharides




Columns: Rezex RCM-Monosaccharide
Aminex HPX-87C
Dimensions: 300 x 7.8 mm
Mobile Phase: Water
Flow Rate: 0.6 mL/min
Detection: ELSD
Temperature: 80 °C
Sample: 1. Melezitose 4. Mannose
2. Maltose 5. Fructose
3. Glucose 6. Ribitol

► Orange Juice




Columns: Rezex RCM-Monosaccharide
Aminex HPX-87C
Dimensions: 300 x 7.8 mm
Mobile Phase: Water
Flow Rate: 0.6 mL/min
Detection: ELSD
Temperature: 80 °C
Sample: 1. Sucrose
2. Glucose
3. Fructose

► Organic Acids in Wine



Columns: Rezex ROA-Organic Acid
Aminex HPX-87H
Dimensions: 300 x 7.8 mm
Mobile Phase: 0.005 N H₂SO₄
Flow Rate: 0.5 mL/min
Detection: UV @ 210 nm
Temperature: 40 °C
Sample: 1. Citric acid 4. Lactic acid
2. Tartaric acid 5. Acetic acid
3. Succinic acid

► Sugar Alcohols





Columns: Rezex RPM-Monosaccharide
Aminex HPX-87P
Dimensions: 300 x 7.8 mm
Mobile Phase: Water
Flow Rate: 0.6 mL/min
Detection: ELSD
Temperature: 80 °C
Sample: 1. Glucose 4. Xylitol
2. Mannitol 5. Erythritol
3. Sorbitol

Aminex and Bio-Rad are registered trademarks of Bio-Rad Laboratories, Inc. Comparative separations may not be representative of all applications.

Phenomenex® Rezex™ vs. Bio-Rad® Aminex®



Lower Backpressure = Longer Column Lifetimes and Faster Run Time Capability

▶ Aliphatic Acids



Columns: Rezex ROA-Organic Acid
Aminex HPX-87H
Dimensions: 300 x 7.8 mm
Mobile Phase: 0.005 N H₂SO₄
Flow Rate: 0.5 mL/min
Detection: UV @ 210 nm
Temperature: 40 °C
Sample: 1. Oxalic acid 4. Succinic acid
2. Citric acid 5. Formic acid
3. Tartaric acid 6. Acetic acid



▶ Saccharides



Columns: Rezex RPM-Monosaccharide
Aminex HPX-87P
Dimensions: 300 x 7.8 mm
Mobile Phase: Water
Flow Rate: 0.6 mL/min
Detection: ELSD
Temperature: 75 °C
Sample: 1. Melezitose 4. Mannose
2. Maltose 5. Fructose
3. Glucose 6. Ribitol



Higher Efficiency = Baseline Separation of Critical Sample Components

▶ Sugars



Columns: Rezex RCM-Monosaccharide
Aminex HPX-87C
Dimensions: 300 x 7.8 mm
Mobile Phase: Water
Flow Rate: 0.6 mL/min
Detection: ELSD
Temperature: 80 °C
Sample: 1. Sucrose 4. Fructose
2. Glucose 5. Sorbitol
3. Galactose

▶ Oligosaccharides




Columns: Rezex RSO-Oligosaccharide
Aminex HPX-42A
Dimensions: 200 x 10.0 mm
Mobile Phase: Water
Flow Rate: 0.3 mL/min
Detection: ELSD
Temperature: 80 °C
Sample: 1. Malto-Oligosaccharides (Dp1-Dp14)
(Dp refers to degree of polymerization)

Aminex and Bio-Rad are registered trademarks of Bio-Rad Laboratories, Inc. Comparative separations may not be representative of all applications.

Select Applications On Rezex™ HPLC Columns

Food and Beverage


▶ Apple Juice



App ID 5505

Column: Rezex RCM-Monosaccharide
Dimensions: 300 x 7.8 mm
Part No.: 00H-0130-K0
Mobile Phase: Water
Flow Rate: 0.6 mL/min
Detection: RI
Temperature: 75 °C
Sample:
 1. Sucrose
 2. Glucose
 3. Fructose


▶ White Wine



App ID 14970

Column: Rezex ROA-Organic Acid
Dimensions: 300 x 7.8 mm
Part No.: 00H-0138-K0
Mobile Phase: 0.005 N H₂SO₄
Flow Rate: 0.5 mL/min
Detection: UV @ 210 nm
Temperature: 40 °C
Sample:
 1. Citric acid
 2. Tartaric acid
 3. Succinic acid
 4. Lactic acid
 5. Acetic acid

▶ Honey




App ID 5503

Column: Rezex RCM-Monosaccharide
Dimensions: 300 x 7.8 mm
Part No.: 00H-0130-K0
Mobile Phase: Water
Flow Rate: 0.6 mL/min
Detection: RI
Temperature: 75 °C
Sample:
 1. Dp 3
 2. Dp 2
 3. Glucose
 4. Fructose

Fermentation and Bioethanol


▶ Fermentation Broth



App ID 14987, 14988, 14889

Column: Rezex RCM-Monosaccharide
Dimensions: 300 x 7.8 mm
Part No.: 00H-0130-K0
Mobile Phase: Water
Flow Rate: 0.5 mL/min
Detection: ELSD
Temperature: 80 °C
Sample:
 1. Sucrose
 2. Glucose
 3. Galactose
 4. Fructose


▶ Bioethanol



App ID 16650

Column: Rezex ROA-Organic Acid
Dimensions: 300 x 7.8 mm
Part No.: 00H-0138-K0
Mobile Phase: 0.005 N Sulfuric Acid
Flow Rate: 0.6 mL/min
Detection: RI @ 40 °C
Temperature: 60 °C
System: Shimadzu™ Prominence™ LC-20A System
Sample:
 1. Dp4+
 2. Dp3
 3. Maltose
 4. Glucose
 5. Lactic Acid
 6. Glycerol
 7. Acetic Acid
 8. Ethanol

▶ Bioethanol




App ID 16648

Column: Rezex ROA-Organic Acid
Dimensions: 150 x 7.8 mm
Part No.: 00F-0138-K0
Mobile Phase: 0.005 N Sulfuric Acid
Flow Rate: 0.6 mL/min
Detection: RI @ 40 °C
Temperature: 60 °C
System: Shimadzu™ Prominence™ LC-20A System
Sample:
 1. Dp4+
 2. Dp3
 3. Maltose
 4. Glucose
 5. Lactic Acid
 6. Glycerol
 7. Acetic Acid
 8. Ethanol

Select Applications On Rezex™ HPLC Columns

Saccharides


► Monosaccharides with Organic Acids



App ID 14969

Column: Rezex RHM-Monosaccharide
Dimensions: 300 x 7.8 mm
Part No.: 00H-0132-KO
Mobile Phase: 0.1 % formic acid
Flow Rate: 0.5 mL/min
Detection: ELSD
Temperature: 75 °C
Sample: 1. Citric acid
 2. Tartaric acid
 3. Glucose
 4. Fructose
 5. Succinic acid


► Saccharides



App ID 5504

Column: Rezex RCM-Monosaccharide
Dimensions: 300 x 7.8 mm
Part No.: 00H-0130-KO
Mobile Phase: Water
Flow Rate: 0.6 mL/min
Detection: RI
Temperature: 85 °C
Sample: 1. Melezitose
 2. Maltose
 3. Glucose
 4. Mannose
 5. Fructose
 6. Ribitol

► Saccharides




App ID 14966

Column: Rezex RAM-Carbohydrate
Dimensions: 300 x 7.8 mm
Part No.: 00H-0131-KO
Mobile Phase: Water
Flow Rate: 0.6 mL/min
Detection: ELSD
Temperature: 75 °C
Sample: 1. Melezitose
 2. Mannose
 3. Maltose
 4. Fructose
 5. Glucose

Carbohydrates


► Carbohydrates



App ID 14974

Column: Rezex RAM-Carbohydrate
Dimensions: 300 x 7.8 mm
Part No.: 00H-0131-KO
Mobile Phase: Water
Flow Rate: 0.6 mL/min
Detection: ELSD
Temperature: 80 °C
Sample: 1. Sucrose
 2. Glucose
 3. Arabinose
 4. Ribose


► Carbohydrates



App ID 5516

Column: Rezex RCU-USP Sugar Alcohols
Dimensions: 250 x 4.0 mm
Part No.: 00G-0130-D0
Mobile Phase: Water
Flow Rate: 0.2 mL/min
Detection: RI
Temperature: 30 °C
Sample: 1. Mannitol
 2. Sorbitol

► Carbohydrates




App ID 5510

Column: Rezex RNM-Carbohydrate
Dimensions: 300 x 7.8 mm
Part No.: 00H-0136-KO
Mobile Phase: Water
Flow Rate: 0.6 mL/min
Detection: RI
Temperature: 75 °C
Sample: 1. Stachyose
 2. Cellobiose
 3. Glucose
 4. Fructose
 5. Arabinose
 6. Ribose
 7. Salicin (1 % solution of each)

Select Applications On Rezex™ HPLC Columns

Carbohydrates / Saccharides (Cont'd)


▶ Carbohydrates



App ID 5506

Column: Rezex RKP-Potassium
Dimensions: 300 x 7.8 mm
Part No.: 00H-3252-KO
Mobile Phase: Water
Flow Rate: 0.4 mL/min
Detection: RI
Temperature: 85 °C
Sample: 1. Maltotriose
 2. Maltose
 3. Glucose


▶ Wood Saccharides



App ID 14967

Column: Rezex RPM-Monosaccharide
Dimensions: 300 x 7.8 mm
Part No.: 00H-0135-KO
Mobile Phase: Water
Flow Rate: 0.6 mL/min
Detection: ELSD
Temperature: 85 °C
Sample: 1. D-Cellobiose
 2. α-D-glucose
 3. D-(+)-Xylose
 4. D-(+)-Galactose
 5. D-(-)-Arabinose

▶ Saccharides and Alcohols




App ID 5508

Column: Rezex RPM-Monosaccharide
Dimensions: 300 x 7.8 mm
Part No.: 00H-0135-KO
Mobile Phase: Water
Flow Rate: 0.6 mL/min
Detection: RI
Temperature: 75 °C
Sample: 1. Stachyose 8. Mannitol
 2. Maltose 9. Salicin
 3. Glucose 10. Xylitol
 4. Xylose 11. Sorbitol
 5. Galactose
 6. Fructose
 7. meso-Erythritol

Alcohols


▶ Alcohols



App ID 5512

Column: Rezex RHM-Monosaccharide
Dimensions: 300 x 7.8 mm
Part No.: 00H-0132-KO
Mobile Phase: Water
Flow Rate: 0.6 mL/min
Detection: RI
Temperature: 60 °C
Sample: 1. Methanol
 2. Ethanol
 3. Isopropanol
 4. n-Propanol
 5. sec-Butanol
 6. n-Butanol


▶ Carboxylic Acids



App ID 5517

Column: Rezex ROA-Organic Acid
Dimensions: 300 x 7.8 mm
Part No.: 00H-0138-KO
Mobile Phase: Water + 0.5 % Trifluoroacetic Acid
Flow Rate: 1.0 mL/min
Detection: RI
Temperature: 40 °C
Sample: 1. Acetylene Carboxylic Acid
 2. Maleic Acid
 3. Succinic Acid
 4. Fumaric Acid

▶ Amino Sugars




App ID 5514

Column: Rezex ROA-Organic Acid
Dimensions: 300 x 7.8 mm
Part No.: 00H-0138-KO
Mobile Phase: 1 % Phosphoric Acid
Flow Rate: 0.6 mL/min
Detection: RI
Temperature: Ambient
Sample: 1. Glucose
 2. N-Acetylglucose
 3. N-Acetylgalactosamine

Select Applications On Rezex™ HPLC Columns


Other Applications

► Food Softeners




Column: Rezex RCM-Monosaccharide
Dimensions: 300 x 7.8 mm
Part No.: 00H-0130-KO
Mobile Phase: Water
Flow Rate: 0.5 mL/min
Detection: RI
Temperature: 60 °C
Sample: 1. Glycerol
 2. Methoxypolyethylene Glycol
 3. Triethylene Glycol
 4. Sorbitol
 5. Urea

► Oligosaccharides



Column: Rezex RSO-Oligosaccharide
Dimensions: 200 x 10 mm
Part No.: 00P-0133-NO
Mobile Phase: Water
Flow Rate: 0.3 mL/min
Detection: RI
Temperature: 75 °C
Sample: Malto-Oligosaccharides as shown

► Oligosaccharides



Column: Rezex RPM-Monosaccharide
Dimensions: 100 x 7.8 mm
Part No.: 00D-0135-KO
Mobile Phase: Water
Flow Rate: 0.6 mL/min
Detection: ELSD (Ambient)
Column Temperature: 80 °C
Sample: 1. Mannitol
 2. Sorbitol

Applications (partial list)


Compound	Rezex Phase	App ID No.
Acetic acid	ROA	5515
	RHM	14962
Acetylene carboxylic acid	ROA	5517
Arabinose	RNM	5510
	RAM	14974
Citric acid	ROA	5515
	RHM	14970
D-Cellobiose	RNM	5510
	RPM	14967
D-(+)-Galactose	RPM	14967
D-(+)-Xylose	RPM	14967
D-(-)-Arabinose	RPM	14967
Erythritol	RPM	14968
Ethanol	ROA	5512
Formic acid	ROA	5515
	RHM	14962
Fructose	RCM	5504, 5505, 14975, 14981, 14977, 14987
	RPM	5508, 14985
	RNM	5510
	RAM	14966
Fumaric acid	ROA	5517
Galactose	RPM	5508
	RCM	14977, 14987
Glucose	RCM	5504, 5505, 14975, 14977, 14979, 14981, 14987, 14983
	RPM	5508, 14985
	RKP	5506
	ROA	5514
	RNM	5510
	RAM	14966, 14974
Glycerol	RCM	5518
Isopropanol	ROA	5512
Lactic acid	RHM	14970, 14971
Maleic acid	ROA	5517
Malto-Oligosaccharides	RSO	5507
Maltose	RCM	5504, 14975
	RPM	5508, 14985
	RKP	5506
	RAM	14966
Maltotriose	RKP	5506
Mannitol	RPM	5508, 14968
	RCU	5516


Compound	Rezex Phase	App ID No.
Mannose	RCM	5504, 14975
	RPM	14985
	RAM	14966
Melezitose	RCM	5504, 14975
	RPM	14985
	RAM	14966
meso-Erythritol	RPM	5508
Methanol	ROA	5512
Methoxypolyethylene Glycol	RCM	5518
n-Butanol	ROA	5512
n-Propanol	ROA	5512
N-Acetylgalactosamine	ROA	5514
N-Acetylglucose	ROA	5514
Oxalic acid	ROA	5515
	RHM	14962
Ribose	RNM	5510
	RAM	14974
Ribitol	RCM	5504, 14975
	RAM	14966
	RHM	14991
	RPM	14985
Salicin	RPM	5508
	RNM	5510
sec-Butanol	ROA	5512
Sorbitol	RPM	5508, 14968
	RCM	5518, 14977
	RCU	5516
Stachyose	RPM	5508
	RNM	5510
Succinic acid	ROA	5515
	RHM	14962
Sucrose	RCM	5505, 14981, 14977, 14987
	RAM	14974
Tartaric acid	RCM	5504
	RHM	14970
Triethylene glycol	RCM	5518
Urea	RCM	5518
Xylitol	RPM	5508, 14968
Xylose	RPM	5508

Call us for assistance in selecting the right column for your needs

* For applications not in this brochure, copies are available through Phenomenex or your local Phenomenex distributor.

Specifications And Operating Recommendations

	RCM Monosaccharide	RSO Oligosaccharide	RNO Oligosaccharide	RNM Carbohydrate	RAM Carbohydrate
Part No.	00H-0130-K0	00P-0133-N0	00P-0137-N0	00H-0136-K0	00H-0131-K0
Ionic Form	Calcium	Silver	Sodium	Sodium	Silver
Standard Dimensions	300 x 7.8 mm	200 x 10 mm	200 x 10 mm	300 x 7.8 mm	300 x 7.8 mm
Matrix	Sulfonated Styrene Divinyl Benzene				
Cross Linking	8 %	4 %	4 %	8 %	8 %
Particle Size	8 µm	12 µm	12 µm	8 µm	8 µm
Min. Efficiency (p/m) based on last peak	35,000	N/A	N/A	30,000	35,000
Typical Pressure (psi @ Max Flow Rate)	350	200	200	350	350
Max. Pressure (psi @ Max Flow Rate)	1,000	300	300	1,000	1,000
Max. Flow Rate (mL/min)	0.6 (see pressure)	0.3	0.3	0.6	0.6
Max. Temperature (°C)	85	85	85	85	85
Typical Mobile Phase	Water	Water	Water	Water	Water
pH Range	Neutral	Neutral	Neutral	Neutral	Neutral
Guard Column Part No.	03B-0130-K0	03R-0133-N0	03R-0137-N0	03B-0136-K0	03B-0131-K0
Cleaning, Regeneration and Storage					
Organic Modifiers (Max)	10 % Methanol, IPA, EtOH, Acetonitrile				
Inorganic Modifiers	5 % CaSO ₄ , Ca(NO ₃) ₂ , CaCl ₂	5 % Silver Nitrate	5 % Sodium Salts	5 % Sodium Salts	2 % Silver Nitrate
Avoid 	>10 % Organic, Acids, Bases, Non-Calcium Salts/ Metal Ions	>5 % Organic, Acids, Bases, Non-Silver Salts/ Metal Ions	>10 % Organic, Acids, Bases, Non-Sodium Salts/ Metal Ions	>10 % Organic, Acids, Bases, Non-Sodium Salts/ Metal Ions	>5 % Organic, Acids, Bases, Non-Silver Salts/ Metal Ions
Cleaning Solvent	100 % Water	100 % Water	100 % Water	100 % Water	100 % Water
Flow Rate(mL/min)	0.6	0.2	0.2	0.4	0.4
Temperature (°C)	85	75	75	75	75
Duration (hrs)	12	12	12	12	12
Regeneration Solvent	0.1 M Ca(NO ₃) ₂	0.1 M AgNO ₃	0.1 M NaNO ₃	0.1 M NaNO ₃	0.1 M AgNO ₃
Flow Rate (mL/min)	0.2	0.1	0.2	0.2	0.2
Temperature (°C)	85	-	85	85	85
Duration (hrs)	4-16	-	4-16	4-16	4-16
Ship/Storage Solvent	Water	Water	Water	Water	Water

	RPM Monosaccharide	RHM Monosaccharide	ROA Organic Acid	RFQ Fast Acid	RCU Sugar Alcohols
Part No.	00H-0135-K0	00H-0132-K0	00H-0138-K0	00D-0223-K0	00G-0130-D0
Ionic Form	Lead	Hydrogen	Hydrogen	Hydrogen	Calcium
Standard Dimensions	300 x 7.8 mm	300 x 7.8 mm	300 x 7.8 mm	100 x 7.8 mm	250 x 4.0 mm
Matrix	Sulfonated Styrene Divinyl Benzene				
Cross Linking	8 %	8 %	8 %	8 %	8 %
Particle Size	8 µm	8 µm	8 µm	8 µm	8 µm
Min. Efficiency (p/m) based on last peak	35,000	35,000	50,000 (Acetic Acid)	30,000	12,000
Typical Pressure (psi @ Max Flow Rate)	350	350	350	350	350
Max. Pressure (psi @ Max Flow Rate)	1,000	1,000	1,000	1,000	1,000
Max. Flow Rate (mL/min)	0.6	0.6	0.6	1.0	0.4
Max. Temperature (°C)	85	85	85	85	85
Typical Mobile Phase	Water	Water	0.005 N H ₂ SO ₄	0.005 N H ₂ SO ₄	Water
pH Range	Neutral	1-8	1-8	1-8	Neutral
Guard Column Part No.	03B-0135-K0	03B-0132-K0	03B-0138-K0	03B-0223-K0	03A-0130-D0
Cleaning, Regeneration and Storage					
Organic Modifiers (Max)	10 % Methanol, IPA, EtOH, Acetonitrile				
Inorganic Modifiers	5 % Lead Nitrate	5 % HNO ₃ , H ₂ PO ₄	5 % HNO ₃ , H ₂ PO ₄	5 % HNO ₃ , H ₂ PO ₄	5 % CaSO ₄ , Ca(NO ₃) ₂ , CaCl ₂
Avoid 	>10 % Organic, Acids, Bases, Non-Lead Salts/ Metal Ions	>10 % Organic, Acids, Bases, Salts, Metal Ions	>10 % Organic, Acids, Bases, Salts, Metal Ions, pH > 3	>10 % Organic, Acids, Bases, Salts, Metal Ions, pH > 3	>10 % Organic, Acids, Bases, Non-Calcium Salts or Metal Ions
Cleaning Solvent	100 % Water	100 % Water	100 % Water	100 % Water	100 % Water
Flow Rate(mL/min)	0.6	0.6	0.5	0.6	0.2
Temperature (°C)	85	85	85	85	85
Duration (hrs)	12	12	12	12	12
Regeneration Solvent	0.1 M Pb(NO ₃) ₂	0.025 M H ₂ SO ₄	0.025 M H ₂ SO ₄	0.025 M H ₂ SO ₄	0.1 M Ca(NO ₃) ₂
Flow Rate (mL/min)	0.2	0.2	0.2	0.2	0.2
Temperature (°C)	85	85	85	85	85
Duration (hrs)	4-16	4-16	4-16	4-16	4-16
Ship/Storage Solvent	Water	Water	0.005 N H ₂ SO ₄	0.005 N H ₂ SO ₄	Water

Ordering Information **REZEX™**

SecurityGuard™ Analytical Cartridges require universal holder Part No.: KJO-4282

Rezex™ Columns						Guards			SecurityGuard™ Cartridges (mm)
Description	Part No.	Cross Linkage	Ionic Form	Size (mm)	Price	Part No.	Size (mm)	Price	4 x 3.0 /10pk
RCM-Monosaccharide	00H-0130-K0	8 %	Calcium	300 x 7.8		03B-0130-K0	50 x 7.8		AJO-4493
RHM-Monosaccharide	00H-0132-K0	8 %	Hydrogen	300 x 7.8		03B-0132-K0	50 x 7.8		AJO-4490
RAM-Carbohydrate	00H-0131-K0	8 %	Silver	300 x 7.8		—	—		AJO-4491
RSO-Oligosaccharide	00P-0133-N0	4 %	Silver	200 x 10.0		03R-0133-N0	60 x 10.0		—
RNO-Oligosaccharide	00P-0137-N0	4 %	Sodium	200 x 10.0		03R-0137-N0	60 x 10.0		—
RPM-Monosaccharide (for USP procedure)	00H-0135-K0	8 %	Lead	300 x 7.8		03B-0135-K0	50 x 7.8		AJO-4492
	00D-0135-K0	8 %	Lead	100 x 7.8		03B-0135-K0	50 x 7.8		AJO-4492
RNM-Carbohydrate	00H-0136-K0	8 %	Sodium	300 x 7.8		03B-0136-K0	50 x 7.8		—
ROA-Organic Acid	00H-0138-K0	8 %	Hydrogen	300 x 7.8		03B-0138-K0	50 x 7.8		AJO-4490
ROA-Organic Acid	00F-0138-K0	8 %	Hydrogen	150 x 7.8		03B-0138-K0	50 x 7.8		AJO-4490
ROA-Organic Acid	00G-0138-E0	8 %	Hydrogen	250 x 4.6		—	—		AJO-4490
RKP-Potassium	00H-3252-K0	8 %	Potassium	300 x 7.8		—	—		—
RFQ-Fast Acid	00D-0223-K0	8 %	Hydrogen	100 x 7.8		03B-0223-K0	50 x 7.8		AJO-4490
RCU-USP Sugar Alcohols	00G-0130-D0	8 %	Calcium	250 x 4.0		03A-0130-D0	30 x 4.0		AJO-4493

for ID: 3.2-8.0 mm



Column Cross Reference Chart

Phenomenex® Rezex™	Bio-Rad® Aminex®	Supelco® SUPELCOGEL™	Waters® Sugar-Pak™
RCM-Monosaccharide	HPX-87C 125-0095	Supelcogel Ca	Sugar-Pak 1
RHM-Monosaccharide	HPX-87H 125-0140	Supelcogel C-610H & H	N/A
RPM-Monosaccharide	HPX-87P 125-0098	Supelcogel Pb	N/A
RNM-Carbohydrate	HPX-87N 125-0143	N/A	N/A
RSO-Oligosaccharide	HPX-42A 125-0097	Supelcogel Ag1 & Ag2	N/A
ROA-Organic Acid	HPX-87H 125-0140	Supelcogel C-610H & H	N/A
RFQ-Fast Acid	Fast Acid 125-0100	N/A	N/A
RKP-Potassium	HPX-87K 125-0142	Supelcogel K	N/A
RCU-USP Sugar Alcohols	Sugar Alcohols 125-0094	N/A	N/A



If you are not completely satisfied with the performance of any Rezex column, as compared to a competing product of the same size and phase, simply return the Rezex column with your comparative data within 45 days for a FULL REFUND.

Phenex™ Syringe Filters

For Sample and Solvent Filtration Prior to Chromatography!

- Less system downtime
- More consistent, reproducible results
- Increased column lifetime



If Phenex™ Syringe Filters do not perform as well or better than your current syringe filter product of similar membrane, diameter and pore size, send in your comparative data within 45 days and keep the Phenex products for FREE!

Syringe filters are non-sterile. Housing is made of medical-grade polypropylene (PP) and offer luer lock inlet/slip outlet connections, unless otherwise indicated.

- 17 mm diameter.
- Glass fiber filters are 28 mm diameter and made of borosilicate. They will remove 90 % of all particles >1.2 µm.
- Housing material is methacrylate butadiene styrene (MBS) polymerisate. Also known as Cyrolite™.
- 26 mm diameter.
- Hydrophobic membrane. Can be made hydrophilic by pre-wetting with IPA.
- 28 mm diameter.
- Additional dimensions and membrane types are available. Please contact your local Phenomenex technical consultant or distributor for availability or assistance.
- Larger quantity purchases at significant savings are available.

Membrane Type/Size	15-17 mm Diameter for 2 – 10 mL sample volumes			25-28 mm Diameter for 10 – 100 mL sample volumes		
	Part No.	Unit	Price	Part No.	Unit	Price
0.45 µm						
Phenex-RC (Regenerated Cellulose)	AF0-2103-12	100/pk		AF0-8103-12 ⁴	100/pk	
	AF0-2103-52	500/pk		AF0-8103-52 ⁴	500/pk	
Phenex-PES ³ (Polyethersulfone)	AF2-5108-12 ¹	100/pk		AF0-8108-12 ⁶	100/pk	
	—	—		AF0-8108-52 ⁶	500/pk	
Phenex-PTFE ⁵ (Polytetrafluoroethylene)	AF0-2102-12	100/pk		AF0-1102-12	100/pk	
	AF0-2102-52	500/pk		AF0-1102-52	500/pk	
Phenex-NY (Nylon)	AF2-5107-12 ¹	100/pk		AF0-1107-12	100/pk	
	AF2-5107-52 ¹	500/pk		AF0-1107-52	500/pk	
0.20 µm						
Phenex-RC (Regenerated Cellulose)	AF0-2203-12	100/pk		AF0-8203-12 ⁴	100/pk	
	AF0-2203-52	500/pk		AF0-8203-52 ⁴	500/pk	
Phenex-PES ³ (Polyethersulfone)	—	—		AF0-8208-12 ⁶	100/pk	
	—	—		AF0-8208-52 ⁶	500/pk	
Phenex-PTFE ⁵ (Polytetrafluoroethylene)	AF0-2202-12	100/pk		AF0-1202-12	100/pk	
	AF0-2202-52	500/pk		AF0-1202-52	500/pk	
Phenex-NY (Nylon)	AF2-5207-12 ¹	100/pk		AF0-1207-12	100/pk	
	AF2-5207-52 ¹	500/pk		AF0-1207-52	500/pk	
1.20 µm						
Phenex-GF ^{2,3} (Glass Fiber)	<i>Prefiltration of heavily contaminated or highly viscous samples. When used in-series preceding a membrane filter, clogging of the membrane filter is prevented and sample clean up is optimized. Outlet connection is luer lock.</i>			AF0-8515-12 ⁶	100/pk	
				AF0-8515-52 ⁶	500/pk	



You Need Constant Temperature, Use A ThermaSphere™ TS-130 Column Heater!

BENEFITS

- Improves chromatographic reproducibility
- Improves ruggedness of separation
- Improves peak efficiency and analyte quantification
- Improves baseline and overall detector performance

ESSENTIAL FOR

- Monosaccharides and sugar alcohol analysis
- Organic acid analysis
- Oligosaccharide analysis
- Rapid Fruit Sugar Analysis

FEATURES

- Compact and lightweight unit
- Transparent compartment cover
- Micro-fan provides rapid thermal equilibration

Specifications

Column Size Accommodated: Fits up to one 30 cm length column, or 25 cm column with guard column. Multiple inlet and outlet slots allow the shortest length of tubing to be used with any length column.

Temperature Range: From 25 to 90 °C in 0.1 °C increments.

Temperature Stability: ±0.1 °C Calibration two-point, electronic, factory set.

Accuracy: 0.5 °C over the entire range.

Power: 12 volt DC universal power supply takes voltage inputs from 95 VAC to 265 VAC, 50/60 Hz. CE approved.

Over-temperature Alarm: Audible with automatic heater shutoff if column temperature exceeds 10 °C of target temperature.

Auto-Off Timer: Count down timer with audible alarm turns off heater, settable to 30 days in days, hours, minutes and seconds.

Injection Counter: Trigger on external switch closure.

ThermaSphere TS-130 Column Heater

Part No.	Description	Price
EHO-7057	ThermaSphere TS-130 HPLC Column Heater 25-90 °C, 95 to 265 VAC, 50/60 Hz	
EHO-7058	Stand for ThermaSphere TS-130 HPLC Column Heater	

1. The ThermaSphere TS-130 is warranted for one year parts and labor. Each unit is individually calibrated and comes with a Certificate of Performance. No adjustment or re-calibration is ever necessary. CE approved, UL and CSA approved power supplies.
2. Please specify Line Cord if other than North America (Australia, Germany, Italy and U.K. are available)



If you are not completely satisfied with the performance of your new TS-130 column heater, simply return the unit within 60 days of purchase for a FULL REFUND.

Trademarks

Phenex, Rezex, ThermaSphere, and SecurityGuard are trademarks of Phenomenex, Inc. Cyrolite is a trademark of CY/RO Industries. Aminex and Bio-Rad are registered trademarks of Bio-Rad Laboratories, Inc. Waters is a registered trademark and Sugar-Pak is a trademark of Waters Corporation. Supelco is a registered trademark and SUPELCOGEL is a trademark of Sigma-Aldrich Biotechnology LP. Transgenomic is a registered trademark and CARBOsep is a trademark of Transgenomic, Inc.

Disclaimer

Comparative separations may not be representative of all applications. Subject to Phenomenex Standard Terms & Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

© 2009 Phenomenex, Inc. All rights reserved.

Australia

t: 02-9428-6444
f: 02-9428-6445
auinfo@phenomenex.com

Austria

t: 01-319-1301
f: 01-319-1300
anfrage@phenomenex.com

Belgium

t: +31 (0)30-2418700
f: +31 (0)30-2383749
beinfo@phenomenex.com

Canada

t: (800) 543-3681
f: (310) 328-7768
info@phenomenex.com

Denmark

t: 4824 8048
f: 4810 6265
dkinfo@phenomenex.com

France

t: 01 30 09 21 10
f: 01 30 09 21 11
franceinfo@phenomenex.com

Germany

t: 06021-58830-0
f: 06021-58830-11
anfrage@phenomenex.com

Ireland

t: 01 247 5405
f: +44 1625-501796
eireinfo@phenomenex.com

Italy

t: 051 6327511
f: 051 6327555
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

Netherlands

t: 030-2418700
f: 030-2383749
nlinfo@phenomenex.com

New Zealand

t: 09-4780951
f: 09-4780952
nzinfo@phenomenex.com

Puerto Rico

t: (800) 541-HPLC
f: (310) 328-7768
info@phenomenex.com

United Kingdom

t: 01625-501367
f: 01625-501796
ukinfo@phenomenex.com

**All other countries:
Corporate Office USA** 

t: (310) 212-0555
f: (310) 328-7768
info@phenomenex.com



www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

411 Madrid Avenue
Torrance, CA 90501-1430
USA