Supplementary protocol



NucleoSpin® 8 PCR Clean-up – centrifuge processing (Rev. 01)

This protocol is only a supplement to the kit's general user manual. Please refer to the kit manual for more detailed information regarding safety instructions, product-specific disclaimers, and especially preparations needed before starting the procedure. The latest version of the user manual is available at www.mn-net.com/usermanuals or can be requested from our technical service (tech-bio@mn-net.com). Material safety data sheets (MSDS) can be downloaded from www.mn-net.com/MSDS.

For the handling of NucleosSpin® 8 PCR Clean-up kits using a centrifuge, additional equipment is necessary:

- A <u>microtiterplate centrifuge</u> which is able to accommodate the NucleoSpin® PCR Clean-up Binding Strips stacked on a MN Square-well Block or Rack of Tube Strips and reaches accelerations of 5,600–6,000 x g (bucket height: 85 mm)
- Starter Set C (REF 740684), containing Column Holders C, Dummy Strips, MN Square-well Blocks, and Rack of Tube Strips (for detailed information, please refer to the Starter Set C manual)
- Round-well Block (REF 740671) and MN Square-well Block (REF740476)
- For transfer of the samples from the Round-well Block to the NucleoSpin® PCR Clean-up Binding Strips, we recommend use of an <u>electronic eight-channel pipetting device</u> with extra long tips capable of holding more than 500 μL.

1 Adjust the volume of reaction mixture

Before starting the purification procedure, add Tris buffer (10 mM, pH 7.0) to adjust the PCR reaction mixture to a final volume of 100 μ L.

Note: Removal of mineral oil is not necessary.

Transfer 100 μ L sample to each well of the Round-well Block (not provided in the kit).

2 Dispense binding buffer to the NucleoSpin® PCR Clean-up Binding Strips

Add 200 μ L of Buffer NT to 100 μ L sample in into the wells of the Round-well Block and pipette up and down several times for mixing.

3 Transfer PCR samples to the NucleoSpin® PCR Clean-up Binding Strips

Insert the NucleoSpin® PCR Clean-up Binding Strips into the Column Holder C. Place Column Holder C on a MN Square-well Block (not supplied with the kit) and transfer all of the sample/ Buffer NT mixture into the wells of the NucleoSpin® PCR Clean-up Binding Strips. Do not moisten the rims while dispensing samples. Moistened rims may cause cross contamination during centrifugation steps.

4 Bind DNA to silica membrane

Place MN Square-well Block and the Column Holder C with the NucleoSpin® PCR Clean-up Binding Strips onto the centrifuge carrier and place it into the robot buckets. Centrifuge at $5,600 \times g$ for 2 min.

Typically, samples will pass through the columns within ≤ 1 min.

5 Wash silica membrane

1st wash

Empty the Block and add 900 μ L of Buffer NT3 to each well of the NucleoSpin® PCR Clean-up Binding Strips and centrifuge again at 5,600 x g for 1–2 min. After centrifugation, discard flow-through collected in the MN Square-well Block.

2nd wash

Repeat wash step with 900 μ L of Buffer NT3. Centrifuge again at 5,600 x g for 1–2 min and discard the flow-through.

6 Dry NucleoSpin® PCR Clean-up Binding Strips

Centrifuge for 5-10 min at $5,600 \times g$ in order to remove residual washing buffer from the silica membrane and for drying the membrane.

7 Elute plasmid DNA

Place the Column Holder C with the NucleoSpin® PCR Clean-up Binding Strips on top of the Rack of Tube Strips (supplied with the kit).

Note: The Rack of Tube Strips has to be filled with tube strips completely for centrifugation.

Dispense **50–150 µL Elution Buffer NE** to each well of the NucleoSpin® PCR Clean-up Binding Strips. Dispense buffer directly onto the membrane. Incubate at room temperature for **1 min**. Centrifuge at **5,600–6,000 x** *g* for **2–3 min**.

Optional: Preheat Elution Buffer NE to 70 °C before dispensing. This will increase recovery for PCR products > 1000 bp.