

	<i>Increase Column Lifetime</i>	<i>Remove Particulates</i>	<i>Remove Proteins</i>	<i>Remove Phospholipids</i>	<i>De-Salt</i>	<i>Solvent-Switching</i>	<i>Extract Target Analytes</i>	<i>Concentrate Sample</i>	
<b>Filtration</b>	●	●							Phenex Filters
<b>Protein Precipitation</b>	●	●	●						Impact
<b>Phospholipid Removal/ Protein Precipitation</b>	●	●	●	●					Phree
<b>QuEChERS</b>	●		●	●					roQ
<b>SLE</b>	●	●	●	●	●	●			Novum
<b>Solid Phase Extraction</b>	●	●	●	●	●	●	●	●	Strata/Strata-X