




NEBNext Ultra Directional RNA Libr Prep for Illumina


Description and Workflow

The NEBNext Ultra RNA Directional Library Prep Kit for Illumina contains enzymes and buffers that are ideally suited for cDNA library preparation for next-generation sequencing.



** Library Validation is performed off-line on 3rd party devices.

Deck Layout



- ① On Deck Thermal Cycler (ODTC)
- ② Cooling module for NEB reagents
- ③ Heated Shaker, cooling position
- ④ Disposable Tips
- ⑤ Reagents
- ⑥ Plates, magnet and primers
- ⑦ Plate stackers, PCR lids
- ⑧ Liquid waste

Recommended Platforms

Platform	Pipetting Channels	96 MPH	Batch size	ODTC *
STAR	8	No	48 samples	Yes
STAR	8	Yes	96 samples	Yes

* On Deck Thermal Cycler

