

New product

July 2017

# VISOCOLOR<sup>®</sup> Powder Pillows free Chlorine

Overview	
Method	Photometric determination of the chlorine content with DPD
Measuring range	0.03–6.00 mg/L Cl <sub>2</sub>
Tests	100 / 1000
Shelf life	5 years from production
REF	936220 / 936220.1
Scope of delivery	100 Powder Pillows in PE bag (936220) 1000 Powder Pillows in box (936220.1)



## Benefits

- Easier dosing of reagents, than with VISOCOLOR<sup>®</sup> ECO. No spoon required
- Very long shelf life of 5 years from production
- Broad measuring range from 0.03 mg/L up to 6.00 mg/L free chlorine
- Opening of Powder Pillows does not require scissors or other tools
- All important information (parameter, LOT, expiry) is printed on every Powder Pillow
- Pre-programmed on compact photometers PF-12<sup>Plus</sup> and PF-3 (version: Drinking Water, Pool)
- **Ready to use in HACH\* photometers**, with pre-programmed HACH\* methods  
Suitable HACH\* methods: 8021, 10245, 10069

## Markets / applications

- Swimming pools
- Drinking water treatment
- Chemical industry

## General information

The addition of chlorine to swimming pools, water reservoirs and water mains is an approved procedure to rid the water of germs. With the correct dose harmful microorganisms are destroyed, many impurities removed and the growth of algae is prevented. However it is imperative that the chlorine content is regularly checked, since excessive chlorine not only impairs the smell and taste of the water but can be hazardous.

One distinguishes between free chlorine and bound chlorine (chloroamines); the sum of both is called total chlorine, which can be detected with VISOCOLOR<sup>®</sup> Powder Pillows total Chlorine (REF 936221 / 936221.1).

## Contact

If you have further questions, please do not hesitate to contact us:

Technical support and customer service	Tel.:	+49 24 21 969-332
	E-Mail:	csc@mn-net.com
For orders (international)	E-Mail:	sales@mn-net.com

\* HACH is a registered trademark of HACH company, Loveland, Colorado. The VISOCOLOR<sup>®</sup> product is neither endorsed, nor sponsored by or affiliated with the HACH Company

